



Features:

- Audible alarm
- Static refrigerated shelves
- Door lock
- Digital display with built in alarms
- Bottom mounted system
- HI/Lo read out
- Remote alarm contacts



Temperature

-18°C to -22°C



Capacity

350 L

570 L



Dimensions

595w x 565d x 1850h

775w x 750d x 1850h

Model	MLF400	MLF600
Gross capacity	350 L	570 L
External dimensions in mm (w / d / h including castors or legs)	595 x 650 x 1850	775 x 750 x 1850
Internal dimensions in mm (w / d / h)	455 x 445 x 1685	635 x 545 x 1685
Packed dimensions in mm (w / d / h including pallet)	640 x 670 x 2010	830 x 790 x 2010
Gross weight -kg (packed weight)	79	102
Net weight - kg	74	94
Energy consumption	795 kWh/y	896 kWh/y
Voltage	240V	240V
Temperature range	-18°C to -22°C	-18°C to -22°C
Ambient temperature range	10°C to 32°C	10°C to 32°C
Refrigerant	R290	R290
Cooling system / defrost	static / manual defrost	static / manual defrost
Exterior cabinet finish	enamelled steel	enamelled steel
Inner liner material	ABS	ABS
Door	solid	solid
Door hinge	right	right
Door handle	pull handle	pull handle
Type of controller	digital	digital
Temperature display	external	external
Data logger	optional	optional
Door ajar alarms	optional	optional
High and low read out	visual / audible	visual / audible
Remote alarm contacts (BMS connectivity for alarms only)	yes	yes
(BMS connectivity for temperatures)	optional	optional
Port entry	optional	optional
Battery back-up for controller	optional	optional
Power requirement	10amp GPO	10amp GPO
Watts	205W	220W
Castors	No adjustable levelling legs	No adjustable levelling legs
Adjustable shelves	no	no
Usable shelf area in mm (w / d)	6 x Large shelf 455 x 405 1 x Small shelf 438 x 202	6 x Large shelf 635 x 505 1 x Small shelf 618 x 302

laftech
Contamination Control Solutions

Call: 1300 306 002 www.laftech.com.au

7 Longford Ct Springvale VIC 3171
P: +(61) 03 9562 3355
info@mediline.com.au