

NLMB1006 – PREMIUM QUALITY PHARMACY & LABORATORY REFRIGERATOR

This amazing premium quality pharmacy & laboratory refrigerator, comes complete with all the standard features plus almost every other feature you could want in a pharmacy & laboratory refrigerator.

Some stand out features include, heated glass doors to prevent condensation in high humidity areas and an easy to use, built-in data logger.

This refrigerator provides an even temperature throughout by its fan-forced system.



** Mediline Australia reserves the right to alter specifications without prior notice.

NLMB1006 – PREMIUM QUALITY PHARMACY & LABORATORY REFRIGERATOR



FEATURES	NLMB1006
Gross Capacity	1006 Litres
Internal Dimensions in mm (w/d/h)	1100 x 684 x 1325
External Dimensions in mm (w/d/h including castors)	1220 x 872 x 1865
Packed Dimensions in mm (w/d/h including pallet)	1335 x 945 x 2000
Gross Weight (packed weight)	237 kg
Net Weight	189 kg
Voltage	240V 50Hz
Temperature Range	+2°C to +8°C
Ambient Temperature Range	+10°C to +32°C
Refrigerant	R600
Cooling System / Defrost	Fan Forced / Automatic
Exterior Cabinet Finish	Enamelled Steel
Inner Liner Material	Enamelled Steel
Door	Heated Glass / Self-Closing
Door Hinge	Left Hand / Right Hand
Door Handle	Yes
Type of Controller	Microprocessor
Temperature Display	External / Digital
Data Logger	Built-in
Door Ajar Alarms	Visual / Audible
High and Low Readout	Visual / Audible
Remote Alarm Contacts	
(BMS Connectivity for Alarms Only)	√
(BMS Connectivity for Temperatures)	Optional
Clip on Data Strip	√
Port Entry	√ (25mm)
Battery Backup for Controller	√ (Lasts approx. 8 Hours)
Amps	2.4
Watts	425
Castors	Yes with brake at rear
Adjustable Shelves	10
Usable Shelf Area in mm (w/d)	510 x 620
Light	Yes with separate switch (light comes on automatically when opening door)
Noise	53 dB(A)



NL1321-03/22

** Mediline Australia reserves the right to alter specifications without prior notice.