

Features:

- Fan forced system
- Voltage free connection
- Door lock
- Door ajar alarms
- Interior LED lighting
- Automatic condensate water disposal
- Fan door switch
- Data logger included
- Battery backup for controller
- Digital display with built in alarms
- Hi / Lo read out
- Self closing door



Temperature

2°C to 8°C

Capacity

- 29 L
- 59 L
- 125 L
- 300 L
- 400 L
- 500 L

Dimensions

- 445w x 458d x 535h
- 450w x 460d x 875h
- 595w x 610d x 840h
- 595w x 611d x 1560h
- 595w x 611d x 1820h
- 700w x 600d x 2020h

Specifications

Model	ML29GP / P	ML59GP / P	ML125GP / P
Gross capacity	29 L	59 L	125 L
External dimensions in mm (w / d / h including castors or legs)	445 x 458 x 535	450 x 460 x 875	595 x 610 x 840
Internal dimensions in mm (w / d / h)	361 x 205 x 400	390 x 390 x 720	505 x 442 x 625
Packed dimensions in mm (w / d / h including pallet)	500 x 575 x 590	660 x 705 x 910	660 x 705 x 910
Gross weight -kg (packed weight)	30	41	50
Net weight - kg	27	38	46
Voltage	220V	220V	220V
Temperature range	2°C to 8°C	2°C to 8°C	2°C to 8°C
Ambient temperature range	10°C to 32°C	10°C to 32°C	10°C to 32°C
Refrigerant	R600	R600	R600
Cooling system / defrost	fan forced	fan forced	fan forced
Exterior cabinet finish	enamelled steel	enamelled steel	enamelled steel
Inner liner material	ABS	ABS	ABS
Door	glass / solid	glass / solid	glass / solid
Door hinge	RH	RH	RH
Door handle	pull handle	pull handle	pull handle
Type of controller	digital	digital	digital
Temperature display	external	external	external
Data logger	included	included	included
Door ajar alarms	yes	yes	yes
High and low read out	visual / audible	visual / audible	visual / audible
Remote alarm contacts (BMS connectivity for alarms only)	yes	yes	yes
(BMS connectivity for temperatures)	optional	optional	optional
Port entry	yes	yes	yes
Battery back-up for controller	yes	yes	yes
Power requirements	10amp GPO	10amp GPO	10amp GPO
Castors	no adjustable levelling legs	no adjustable levelling legs	no adjustable levelling legs
Adjustable shelves	yes	yes	yes
Usable shelf area in mm (w / d)	2 x 376 x 190	2 x 368 x 271 1 x 376 x 190	2 x 488 x 411 2 x 488 x 276
Light	LED	LED	LED
Noise - dB(A)	45	45	45

Specifications

Model	ML300GP / P	ML400GP / P	ML500GP / P
Gross capacity	300 L	400 L	500 L
External dimensions in mm (w / d / h including castors or legs)	595 x 611 x 1560	595 x 611 x 1820	700 x 600 x 2020
Internal dimensions in mm (w / d / h)	505 x 442 x 1345	505 x 422 x 1605	610 x 477 x 1735
Packed dimensions in mm (w / d / h including pallet)	660 x 705 x 1630	660 x 705 x 1890	760 x 735 x 2100
Gross weight -kg (packed weight)	77	87	132
Net weight - kg	73	82	122
Voltage	220V	220V	220V
Temperature range	2°C to 8°C	2°C to 8°C	2°C to 8°C
Ambient temperature range	10°C to 32°C	10°C to 32°C	10°C to 32°C
Refrigerant	R600	R600	R600
Cooling system / defrost	fan forced	fan forced	fan forced
Exterior cabinet finish	enamelled steel	enamelled steel	enamelled steel
Inner liner material	ABS	ABS	ABS
Door	glass / solid	glass / solid	glass / solid
Door hinge	RH	RH	RH
Door handle	pull handle	pull handle	pull handle
Type of controller	digital	digital	digital
Temperature display	external	external	external
Data logger	included	included	included
Door ajar alarms	yes	yes	yes
High and low read out	visual / audible	visual / audible	visual / audible
Remote alarm contacts (BMS connectivity for alarms only)	yes	yes	yes
(BMS connectivity for temperatures)	optional	optional	optional
Port entry	yes	yes	yes
Battery back-up for controller	yes	yes	yes
Power requirements	10amp GPO	10amp GPO	10amp GPO
Castors	no adjustable levelling legs	no adjustable levelling legs	no adjustable levelling legs
Adjustable shelves	yes	yes	yes
Usable shelf area in mm (w / d)	6 x 488 x 411 1 x 488 x 276	6 x 488 x 411 1 x 488 x 276	6 x 593 x 470 1 x 593 x 300
Light	LED	LED	LED
Noise - dB(A)	45	45	45

laftech
Contamination Control Solutions

Call: 1300 306 002 www.laftech.com.au

