



MLF115 Spark Safe Medical & Laboratory Freezer provides the reliability that you need for your valuable items.

The static refrigerated shelves help to keep an even temperature throughout. Equipped with door lock and door ajar alarm, the MLF115 gives you that security that you need.

Features:

- Door lock
- Power failure battery back up for controller
- Door ajar alarms
- Voltage free connection
- Digital display with built in alarms
- Hi / Lo read out



Temperature

-18°C to -22°C



Capacity

115 L



Dimensions

595w x 611d x 840h

Model	MLF115
Gross capacity	115 L
External dimensions in mm (w / d / h including castors or legs)	595 x 611 x 840
Internal dimensions in mm (w / d / h)	455 x 405 x 575
Packed dimensions in mm (w / d / h including pallet)	660 x 705 x 910
Gross weight - kg (packed weight)	55
Net weight- kg	60
Voltage	220V
Temperature range	-18°C to -22°C
Ambient temperature range	10°C to 32°C
Refrigerant	R404a
Cooling system / defrost	static / manual defrost
Exterior cabinet finish	enamelled steel
Inner liner material	ABS
Door	solid
Door hinge	RH
Door handle	pull handle
Type of controller	digital
Temperature display	external
Data logger	optional
Door ajar alarms	visual / audible
High and low read out	visual / audible
Remote alarm contacts (BMS connectivity for alarms only)	yes
(BMS connectivity for temperatures)	optional
Port entry	no
Battery back-up for controller	yes
Power requirement	10amp GPO
Watts	420w
Castors	no adjustable levelling legs
Adjustable shelves	yes
Usable shelf area in mm (w / d)	3 x 438 x 400

laftech
Contamination Control Solutions

Call: 1300 306 002 www.laftech.com.au