



**MEDI LINE**  
AUSTRALIA

## MLF90 -40 Freezer



The MLF -40 °C Freezer provides the reliability that you need for your valuable items.

The static refrigerated shelves help to keep an even temperature throughout. Equipped with door lock and door ajar alarm, the MLF gives you the security that you need.

### Features:

- Door lock
- Power failure battery back-up for controller
- Door ajar alarms
- Voltage free connection
- Digital display with built in alarms
- Hi/Lo read out
- Port entry
- Built in data logger



### Temperature



- 10°C to - 40°C



### Capacity

90 L



### Dimensions

700w x 669d x 854h

Model	MLF90 -40
Gross capacity	90 L
External dimensions in mm (w / d / h including castors or legs)	700 x 669 x 854
Internal dimensions in mm (w / d / h )	500 x 485 x 532
Packed dimensions in mm (w / d / h including pallet)	765 x 735 x 935
Gross weight -kg (packed weight)	60
Net weight - kg	54
Energy consumption	1.69(kWh/24h)/(AT 25°C )
Voltage	240V
Temperature range	-10°C to -40°C
Ambient temperature range	16°C to -32°C
Refrigerant	R290
Cooling system / defrost	static
Exterior cabinet finish	enamelled steel
Inner liner material	ABS
Door	solid
Door hinge	RH
Door handle	lockable
Type of controller	digital
Temperature display	external
Data logger	yes
Door ajar alarms	yes
High and low read out	visual / audible
Remote alarm contacts (BMS connectivity for alarms only)	yes
(BMS connectivity for temperatures)	optional
Port entry	yes
Battery back-up for controller	yes
Power requirement	10amp GPO
Watts	131W
Castors	2 castors + 2 leveling feet
Adjustable shelves	no - fixed drawers x 3
Usable drawer area in mm (w / d)	1pc (494 x 260 x 156.5) 2pc (494 x 418.5 x 165)

**laftech**  
Contamination Control Solutions

Call: 1300 306 002 [www.laftech.com.au](http://www.laftech.com.au)