



R1500CX

Ventilated refrigeration with uniform temperature, hermetic compressor with air condensation, quiet operation, HFC refrigerant (CFC & HCFC free for environmental safe-guard). Automatic defrosting with automatic evaporation of the condensed water. Tropicalised execution to withstand ambient temperatures up to 43 Deg C



R1500V
Pictured Without Stainless Exterior

- ✓ Choice of glass (C) or solid doors (C)
- ✓ HTS 5 inch LCD Digital controller
- ✓ Fully alarmed
- ✓ Remote Alarm Contacts
- ✓ Adjustable Shelving
- ✓ Lockable
- ✓ Economical to run
- ✓ Reliable quality European built
- ✓ Access port
- ✓ Single body structure with AISI 304 stainless steel with shiny finish and rounded angles for maximum hygiene and cleanliness
- ✓ High density insulation 60mm thickness
- ✓ Reversible doors
- ✓ Magnetic door gaskets and automatic closing



Modern HTS Controller

- ✓ 5 Inch LCD Touch Screen (will operate with gloved hands)
- ✓ MicroSD Data Logging
- ✓ USB port for downloading thermoregulation data
- ✓ Password parameter access
- ✓ Battery back up
- ✓ Energy saving Eco/night mode option

Specifications



APPARECCHI
SCIENTIFICI

Description

Model	KLAB R1500VX/R1500CX HTS
Power	230V 50Hz
Temperature Range	0 - 10 Deg C
Working Temperature	4 Deg C
Capacity	1500L
Exterior Dimensions (WxDxH)	143 x 85 x 203 cm (V) 143 x 80 x 203 cm (C)
Net Weight	220kg (V) 240kg (C)
Main Construction	Stainless Steel Exterior & Interior with insulated Glass Door or solid door
Controller Type	HTS 5"
Number of Doors	2
Number of Shelves	3 X 2 = 6
Access Port	1
Power Consumption (Average)	0.15 Kw
Castors	Y
Lock	Y
Auto Defrost	Y
Remote Alarm Contacts	Y
Certifications (Equipment certified in accordance with)	CEI 66/5 UNI EN 61010-1
	ISO 9001:2008

X denotes AISI 304 Stainless steel V denotes Glass Doors
C denotes solid door

LAF Technologies Pty Ltd

Head Office: 12 Royan Place, Bayswater North, VIC 3153 Australia



Accredited Laboratory

laftech

Call: 1300 306 002 www.laftech.com.au