

# 60L - Corrosive Substance Storage Cabinet

## SC0608

### Standard Features

- ✓ Double walled sheet steel construction with a 40mm thermic air barrier
- ✓ High capacity liquid tight sump to contain spills - 150mm deep
- ✓ Self-closing, close fitting doors held shut by catches at 2 points. The doors do not open inwards and are capable of being opened from within the cabinet  
Built in flash arrestors and vent openings (for mechanical ventilation systems)
- ✓ Strong fully adjustable shelving perforated for free air movement, 100mm adjustment increments
- ✓ Sequential door closing system
- ✓ Stainless steel pin with continuous hinging for strength and reliability
- ✓ State of the art powder coat finish for high durability and performance
- ✓ PVC spill trays on shelves for corrosion protection  
Distinct safety signs and directions that comply with AS 1216 and AS 1319
- ✓ Designed and manufactured in full conformance to AS 3780:2023 - The storage and handling of corrosive substances



### SC0608 Specification

Dimensions (mm)			
	H	W	D
External	1065	500	450
Internal	818	418	370

<b>Shelf Levels</b>	2
<b>Max Capacity</b>	60L
<b>Sump Capacity</b>	23L
<b>Weight (kg)</b>	41.0 kg

