



Tabletop Equipments

Hotplate Stirrer
Vortex Mixer
Overhead Stirrer
Homogenizer

including

Pop-Off Cup Head
3" Rubber Platform
3,300 rpm

JSVX-7R Vortex Mixer

Pop-Off Cup Head



3" Rubber Platform



Interchangeable

Supplied with Interchangeable attachments

- 3 inch Rubber Platform Head
- Pop-Off Cup Head
- Easy Remove Applicator



Functional Design

Suction-Cup type foot provides strong and stable attachment on the table surface and reduces vibration



Aluminum Casting Body

Rugged aluminum casting body with epoxy powder coating
Durable and rust free engineering plastic base



Three Way Switch

Multi-Function Switch


- Continuous ON
- Touch ON
- OFF



Brushless Motor

Brushless shade motor for strong power and long shelf life

#JSVX-7R Vortex Mixer

# Technical Data 	JSVX-7R
Overall Dimensions	W130 x D190 x H120 mm (main body only)
Speed Range Movement	0 to 3,300 rpm Orbital, 4 mm Stroke
Control Movement Selection	Stepless Variable Speed Control 3-Way Switch [Continuous ON / Touch ON / OFF]
Motor Motor Capacity Max. Load	Brushless Shade Motor 30 Watt 0.5 kgs
Safety	Over Current Cut-Off
Materials Body Platform	Aluminum casting with epoxy powder coating Engineering Plastic, Rubber
Electrical Requirement	230VAC 50Hz or 230VAC 60Hz, 1-Phase

PHA-01V 3" Rubber Platform Head



PHA-02V Pop-Off Cup Head



PHA-03V Remove Applicator



Ordering Information

Art. No.	Descriptions
JSVX-7R1	Vortex Mixer 230 VAC 50 Hz
JSVX-7R2	Vortex Mixer 230 VAC 60 Hz
	including 1 ea x 3 inch Rubber Platform Head 1 ea x Pop-Off Cut Head 1 ea x Remove Applicator
Optional Accessories	
PHA-01V	3 inch Rubber Platform Head
PHA-02V	Pop-Off Cup Head
PHA-03V	Remove Applicator