


# The New Bugbox Mx

World's Smallest Anaerobic/Microaerophilic Workstation



N <sub>2</sub>	
CO <sub>2</sub>	0.1 - 30.0 %
O <sub>2</sub>	0.0 - 20.0 %



Width	790 mm
Depth	679 mm
Height	872 mm



①	Optical oxygen sensor
②	Removable front screen
③	SD card Data logging
④	Hypoxic cycling
⑤	Single Plate Entry System (SPES)
⑥	Internal power socket
⑦	Optional internal/external HEPA filter