


SafeDEVELOP™ DFO and Ninhydrin Fingerprint Development Chamber

SD-34S-A - 36"W, 115VAC, 60Hz



DIMENSIONS & WEIGHTS

- Internal Height:** 19.6" | 495 mm
- External Width:** 36.25" | 921 mm
- External Depth:** 25.5" | 648 mm
- External Height:** 29.25" | 743 mm
- Shipping Width:** 40" | 1016 mm
- Shipping Depth:** 48" | 1219 mm
- Net Weight:** 200 lbs | 91 kg
- Shipping Weight:** 245 lbs | 111 kg

CONSTRUCTION

- Cabinet Style:** Desktop
- Construction:** Epoxy-coated Steel

PROTECTION & COMPLIANCE

- Environmental:** ISO 14001 : 2015 ; ENERGY STAR® Partner
- Quality:** ISO 9001 : 2015

OPTIONS AND ACCESSORIES

- Base Stand, Mobile, with Casters:** CART-30
- Vent-Box Filtration*:** VB-24

* Factory installed; specify when ordering

WARRANTY

Air Science Legacy Limited™ Warranty

ELECTRICAL

- Electrical Voltage:** 120V, 60Hz
- NEMA Plug:** 5-20P

CONTROLS & MONITORING

- Control:** Electronic#Programmable Heat and Humidity Control
- Monitoring Options:** // Standard



Air Science® • 120 6th Street • Fort Myers, FL 33907
T/239.489.0024 | 800-306-0656

