


Fume CHAMBER™

BLC-EX-A - 110V, 60Hz



DIMENSIONS & WEIGHTS

External Width: 9.5" | 241 mm

External Depth: 9.5" | 241 mm

External Height: 9.5" | 241 mm

Shipping Width: 16" | 406 mm

Shipping Depth: 16" | 406 mm

Net Weight: 17 lbs | 7 kg

Shipping Weight: 19 lbs | 9 kg

CONSTRUCTION

Construction: Epoxy-coated Steel

PROTECTION & COMPLIANCE

Environmental: ISO 14001 : 2015 ; ENERGY STAR® Partner

Quality: ISO 9001 : 2015

OPTIONS AND ACCESSORIES

Portable Humidifier and Hose: PHUM-BLC

* Factory installed; specify when ordering

WARRANTY

Air Science Legacy Limited™ Warranty

ELECTRICAL

Electrical Voltage: 120V, 60Hz



Air Science® • 120 6th Street • Fort Myers, FL 33907

T/239.489.0024 | 800-306-0656

©2025 Air Science® USA LLC. All Rights Reserved. 03-28-2025 | Specifcations subject to change without notice.