


Fume CHAMBER™

BLC-6.5-A - 110V, 60Hz



DIMENSIONS & WEIGHTS

Internal Height: 31" | 787 mm

External Width: 31" | 787 mm

External Depth: 13" | 330 mm

External Height: 34" | 864 mm

Shipping Width: 36" | 914 mm

Shipping Depth: 18" | 457 mm

Net Weight: 43 lbs | 19 kg

Shipping Weight: 45 lbs | 20 kg

CONSTRUCTION

Construction: Epoxy-coated Steel

PROTECTION & COMPLIANCE

Environmental: ISO 14001 : 2015 ; ENERGY STAR® Partner

Quality: ISO 9001 : 2015

OPTIONS AND ACCESSORIES

Portable Humidifier and Hose: PHUM-BLC

* Factory installed; specify when ordering

WARRANTY

Air Science Legacy Limited™ Warranty

ELECTRICAL

Electrical Voltage: 120V, 60Hz



Air Science® • 120 6th Street • Fort Myers, FL 33907

T/239.489.0024 | 800-306-0656

©2025 Air Science® USA LLC. All Rights Reserved. 03-28-2025 | Specifcations subject to change without notice.