



- NOTES:
- 1.0 CUSTOMER:
  - 2.0 SUMIT DRAWINGS FOR APPROVAL.
  - 3.0 STANDARD FEATURES:
  - 3.1 CONSTRUCTED OF A WELDED FRAME WITH POWDER COATED FINISHING.
  - 3.2 HEAVY DUTY CASTORS.
  - 3.3 LOCKING STORAGE DRAWERS.
  - 3.4 WORK SURFACES, EPOXY RESIN.
  - 3.5 KRAFT PAPERS DISPENSERS.
  - 3.6 ONE FAN FILTERS UNIT.



Designed by CFG	Checked by	Approved by	Date 3/31/2014	Date
AIR SCIENCE USA		EVIDENCE BENCH 72X39X36		
EvidenceBenchR72.Assem			Edition	Sheet 1 / 1