

# NLMF SERIES LABORATORY FREEZER

**NULINE**  
REFRIGERATION

## FEATURES

- Top Mounted System
- Heated Glass and Heated Cabinet Facia to prevent condensation up to 80% R.H.
- Interior Lighting LED
- Automatic Water Disposal
- Digital Display with Built-in Alarms
- Quick temperature recovery system
- ARTG Number 269957



-18°C to -22°C



700 L  
1400 L



770w x 830d x 1990 h  
1480w x 830d x 1990h

**laftech**

Controlled Environment Solutions

Call: 1300 306 002 [www.laftech.com.au](http://www.laftech.com.au)



# NLMF SERIES LABORATORY FREEZER



| Model                                                         | NLMF700/1 G                                 | NLMF700/1 P | NLMF1400/2 G      | NLMF1400/2 P |
|---------------------------------------------------------------|---------------------------------------------|-------------|-------------------|--------------|
| Gross Capacity                                                | 700 L                                       |             | 1400 L            |              |
| External Dimensions in mm<br>(w/d/h including castors)        | 770 x 830 x 1990                            |             | 1480 x 830 x 1990 |              |
| Internal Dimensions in mm<br>(w/d/h)                          | 620 x 630 x 1405                            |             | 1330 x 630 x 1405 |              |
| Packed Dimensions in mm<br>(w/d/h including pallet)           | 840 x 930 x 2140                            |             | 1550 x 930 x 2140 |              |
| Gross Weight - kg (packed weight)                             | 107                                         |             | 137               |              |
| Net Weight - kg                                               | 95                                          |             | 125               |              |
| Voltage                                                       | 240V                                        |             |                   |              |
| Temperature Range                                             | -18°C to -22°C                              |             |                   |              |
| Ambient Temperature Range                                     | 10°C to 32°C                                |             |                   |              |
| Refrigerant                                                   | R404a                                       |             |                   |              |
| Cooling System / Defrost                                      | Fan Forced / Automatic                      |             |                   |              |
| Exterior Cabinet Finish                                       | White Colourbond                            |             |                   |              |
| Inner Liner Material                                          | White Colourbond (optional Stainless Steel) |             |                   |              |
| Door                                                          | Heated Glass                                | Solid       | Heated Glass      | Solid        |
| Door Hinge                                                    | RH                                          |             | LH / RH           |              |
| Door Handle                                                   |                                             |             | Lockable          |              |
| Type of Controller                                            |                                             |             | Digital           |              |
| Temperature Display                                           |                                             |             | External          |              |
| Data Logger                                                   |                                             |             | Optional          |              |
| Door Ajar Alarms                                              |                                             |             | Optional          |              |
| High and Low Readout                                          |                                             |             | Visual / Audible  |              |
| BMS Connectivity<br>(BMS Connectivity for Remote Alarms only) |                                             |             | Optional          |              |
| (BMS Connectivity for 4-20mA Transmitter)                     |                                             |             | Optional          |              |
| Port Entry                                                    |                                             |             | Optional          |              |
| Battery Back-up for Controller                                |                                             |             | N/A               |              |
| Power Requirement                                             |                                             |             | 10amp GPO         |              |
| Castors                                                       |                                             |             | Yes               |              |
| Adjustable Shelves                                            | 5                                           |             | 10                |              |
| Usable Shelf Area in mm (w/d)                                 | 600 x 610                                   |             | 655 x 610         |              |
| Light                                                         | LED                                         |             |                   |              |

Nuline Refrigeration reserves the right to alter the specifications without prior notice.