

NLMF SERIES LABORATORY FREEZER



FEATURES

- Top Mounted System
- Heated Glass and Heated Cabinet Facia to prevent condensation up to 80% R.H.
- Interior Lighting LED
- Automatic Water Disposal
- Digital Display with Built-in Alarms
- Quick temperature recovery system
- TGA Certified



-18°C to -22°C



700 L
1400 L



770w x 830d x 2050h
1480w x 830d x 2050h



Australia

1800 358 101
sales.au@westlab.com
www.westlab.com.au

Canada

1 877 822 1455
sales.ca@westlab.com
www.westlab.com

New Zealand

+64 9 553 3677
sales.nz@westlab.com
www.westlab.co.nz

United States

646 568 5391
sales.us@westlab.com
www.westlab.com

NLMF SERIES LABORATORY FREEZER



Model	NLMF700/1	NLMF1400/2
Gross Capacity	700 L	1400 L
External Dimensions in mm (w/d/h including castors)	770 x 830 x 2050	1480 x 830 x 2050
Internal Dimensions in mm (w/d/h)	620 x 630 x 1465	1330 x 630 x 1465
Packed Dimensions in mm (w/d/h including pallet)	840 x 930 x 2250	1550 x 930 x 2250
Gross Weight - kg (packed weight)	112	137
Net Weight - kg	95	125
Voltage	240V	
Temperature Range	-18°C to -22°C	
Ambient Temperature Range	10°C to 32°C	
Refrigerant	R404a	
Cooling System / Defrost	Fan Forced / Automatic	
Exterior Cabinet Finish	White Colourbond	
Inner Liner Material	White Colourbond (optional Stainless Steel)	
Door	Solid / Heated Glass	
Door Hinge	RH	LH / RH
Door Handle	Lockable	
Type of Controller	Digital	
Temperature Display	External	
Data Logger	Optional	
Door Ajar Alarms	Optional	
High and Low Readout	Visual / Audible	
BMS Connectivity (BMS Connectivity for Remote Alarms only)	Optional	
(BMS Connectivity for 4-20mA Transmitter)	Optional	
Port Entry	Optional	
Battery Back-up for Controller	N/A	
Power Requirement	10amp GPO	
Castors	Yes	
Adjustable Shelves	5	10
Usable Shelf Area in mm(w/d)	600 x 610	655 x 610
Light	LED	

Nuline Refrigeration reserves the right to alter the specifications without prior notice.

NLMF SERIES LABORATORY FREEZER



Since 1985, Nuline Refrigeration has been manufacturing refrigerators and freezer cabinets for the Australian market. Our specialised range of biomedical cabinets are used extensively across the hospitals, medical centres, pharmacies and laboratories.

Focusing on customer needs, we manufacture, import and customise the biomedical cabinets from a range of selected global suppliers. Our factory is based in Melbourne, Australia and we are proudly 100% Australian owned business.

Each of our manufactured cabinets are custom built based on customer orders by our craftsmen right here in Melbourne.

We have extensive capability to design, engineer and build a range of refrigerated cabinets for you. Our standard product range includes, fridges, freezers (up to -40°C), blanket warmer, mortuary refrigerators and refrigerated incubators.

