


EDU™ Ductless Classroom Demonstration Fume Hoods

P10XL-CART-RC-P-A - Mobile fume hood on casters, 34"W, remote controls, 115VAC, 60Hz



DIMENSIONS & WEIGHTS

- Internal Height:** 38" | 965 mm
- External Width:** 34" | 864 mm
- External Depth:** 27.375" | 695 mm
- External Height:** 87" | 2210 mm
- Shipping Width:** 48" | 1219 mm
- Shipping Depth:** 40" | 1016 mm
- Net Weight:** 234 lbs | 107 kg
- Shipping Weight:** 284 lbs | 129 kg

CONSTRUCTION

Construction: Epoxy-coated Steel

PROTECTION & COMPLIANCE

- OSHA:** OSHA Standard -29 CRF, Safety and Health Regulations for General Industry, 1910.1450: Occupational exposure to hazardous chemicals in laboratories. Part B, definition, laboratory type hood. This product may assist you with compliance or as part of your chemical hygiene plan. Please consult your Safety Officer and/or Industrial Hygienist.
- Education:** CLEAPPS Instruction Approved (EDU) based on ASTM-100 filter (UK)
- Environmental:** ISO 14001 : 2015; ENERGY STAR® Partner
- Quality:** ISO 9001 : 2015

OPTIONS AND ACCESSORIES

- Docking Station*:** DOCK-S
- Dual Access*:** DSASH-P10XL
- Duplex Electrical Outlet*:** AS-GFI
- FSA/Autocal Control Panel*:** ADV-P
- Filter Saturation Alarm*:** FSA
- Monitair Controller*:** MON-P
- Polypropylene Construction*:** P10XL-PP-CART-RC
- SafeSwitch HEPA Filter Shutter System:** ASTM-030-SS
- Spill Tray (Stainless):** TRAY-P10XL-SS
- Utilities Package*:** WATER-EDU

* Factory installed; specify when ordering

AIRFLOW & FILTRATION

- Airflow Direction:** Upflow
- Main Filter:** 22

WARRANTY

Air Science Legacy Limited™ Warranty

ELECTRICAL

Electrical Voltage: 120V, 60Hz

NEMA Plug: 5-15P

CONTROLS & MONITORING

Monitoring Options: // Standard



Air Science® • 120 6th Street • Fort Myers, FL 33907

T/239.489.0024 | 800-306-0656

©2025 Air Science® USA LLC. All Rights Reserved. 07-03-2025 | Specifcations subject to change without notice.