

NLMF SERIES LABORATORY FREEZER

NULINE

REFRIGERATION

Est. 1985

FEATURES

- Top Mounted System
- Heated Glass and Heated Cabinet Facia to prevent condensation up to 80% R.H.
- Interior Lighting LED
- Automatic Water Disposal
- Digital Display with Built-in Alarms
- Quick temperature recovery system
- ARTG Number 269957



-18°C to -22°C



700 L
1400 L



770 x 830 x 1990
1480 x 830 x 1990



NLMF SERIES LABORATORY FREEZER



Model	NLMF700/1 G	NLMF700/1 P	NLMF1400/2 G	NLMF1400/2 P
Gross Capacity	700 L		1400 L	
External Dimensions in mm (w/d/h including castors)	770 x 830 x 1990		1480 x 830 x 1990	
Internal Dimensions in mm (w/d/h)	620 x 630 x 1405		1330 x 630 x 1405	
Packed Dimensions in mm (w/d/h including pallet)	840 x 930 x 2140		1550 x 930 x 2140	
Gross Weight - kg (packed weight)	107		137	
Net Weight - kg	95		125	
Voltage	240V			
Temperature Range	-18°C to -22°C			
Ambient Temperature Range	10°C to 32°C			
Refrigerant	R448A			
Cooling System / Defrost	Fan Forced / Automatic			
Exterior Cabinet Finish	White Colourbond			
Inner Liner Material	White Colourbond (optional Stainless Steel)			
Door	Heated Glass	Solid	Heated Glass	Solid
Door Hinge	RH		LH / RH	
Door Handle			Lockable	
Type of Controller			Digital	
Temperature Display			External	
Data Logger			Optional	
Door Ajar Alarms			Optional	
High and Low Readout			Visual / Audible	
BMS Connectivity (BMS Connectivity for Remote Alarms only)			Optional	
(BMS Connectivity for 4-20mA Transmitter)			Optional	
Port Entry			Optional	
Battery Back-up for Controller			N/A	
Power Requirement			10amp GPO	
Castors			Yes	
Adjustable Shelves	5		10	
Usable Shelf Area in mm (w/d)	600 x 610		655 x 610	
Light	LED			

Nuline Refrigeration reserves the right to alter the specifications without prior notice.

AUSTRALIAN MADE

NULINE

REFRIGERATION

EST. 1985

Since 1985, Nuline Refrigeration has been manufacturing refrigerators and freezer cabinets for the Australian market. Our specialised range of biomedical cabinets are used extensively across the hospitals, medical centres, pharmacies and laboratories.

Focusing on customer needs, we manufacture, import and customise the biomedical cabinets from a range of selected global suppliers. Our factory is based in Melbourne, Australia and we are proudly 100% Australian owned business.

Each of our manufactured cabinets are custom built based on customer orders by our craftsmen right here in Melbourne.

We have extensive capability to design, engineer and build a range of refrigerated cabinets for you. Our standard product range includes, fridges, freezers (up to -40°C), blanket warmer, mortuary refrigerators and refrigerated incubators.

