

NLMS SPARK SAFE LABORATORY REFRIGERATOR

NULINE

REFRIGERATION

Est. 1985

FEATURES

- Top Mounted System
- E-glass and heated cabinet facia to prevent condensation up to 80% R.H.
- Glass or Solid Door
- Auto Water Disposal
- Hot Gas Defrost
- Castors
- Door Lock
- Digital Display with Built-in Alarms
- Hi/Lo Readout
- ARTG Number 278162



	2°C to 8°C
	400 L 700 L 1000 L 1614 L
	600 x 620 x 1990 870 x 700 x 1990 1200 x 700 x 1990 1880 x 700 x 1990



NLMS SPARK SAFE LABORATORY REFRIGERATOR



Model	NLMS400/1	NLMS700/1	NLMS1000/2	NLMS1614/3
Gross Capacity	400 L	700 L	1000 L	1614 L
External Dimensions in mm (w/d/h including castors)	600 x 620 x 1990	870 x 700 x 1990	1200 x 700 x 1990	1880 x 700 x 1990
Internal Dimensions in mm (w/d/h)	500 x 470 x 1490	770 x 550 x 1445	1100 x 550 x 1445	1780 x 550 x 1410
Packed Dimensions in mm (w/d/h including pallet)	700 x 720 x 2140	970 x 800 x 2140	1300 x 800 x 2140	1980 x 800 x 2140
Gross Weight - kg (packed weight)	97	107	137	187
Net Weight - kg	85	95	125	175
Voltage	240V			
Temperature Range	2°C to 8°C			
Ambient Temperature Range	10°C to 32°C			
Refrigerant	R448a			
Cooling System / Defrost	Static			
Exterior Cabinet Finish	White Colourbond			
Inner Liner Material	White Colourbond			
Door	Solid / Glass E-Glass			
Door Hinge	RH		LH / RH	LH / RH / RH
Door Handle	Lockable			
Type of Controller	Digital			
Temperature Display	External			
Data Logger	Optional			
Door Ajar Alarms	Optional			
High and Low Readout	Visual / Audible			
BMS Connectivity (BMS Connectivity for Remote Alarms only)	Optional			
(BMS Connectivity for 4-20mA Transmitter)	Optional			
Port Entry	Optional			
Battery Back-up for Controller	N/A			
Power Requirement	10amp GPO			
Castors	Yes			
Adjustable Shelves	5			15
Usable Shelf Area in mm (w/d)	480 x 450	750 x 530	1080 x 530	565 x 530

Nuline Refrigeration reserves the right to alter the specifications without prior notice.

AUSTRALIAN MADE

NULINE

REFRIGERATION

EST. 1985

Since 1985, Nuline Refrigeration has been manufacturing refrigerators and freezer cabinets for the Australian market. Our specialised range of biomedical cabinets are used extensively across the hospitals, medical centres, pharmacies and laboratories.

Focusing on customer needs, we manufacture, import and customise the biomedical cabinets from a range of selected global suppliers. Our factory is based in Melbourne, Australia and we are proudly 100% Australian owned business.

Each of our manufactured cabinets are custom built based on customer orders by our craftsmen right here in Melbourne.

We have extensive capability to design, engineer and build a range of refrigerated cabinets for you. Our standard product range includes, fridges, freezers (up to -40°C), blanket warmer, mortuary refrigerators and refrigerated incubators.

